

## New England Biolabs Product Specification

*Product Name:* pCLIPf Vector  
*Catalog #:* N9215S  
*Concentration:* 500 µg/ml  
*Shelf Life:* 24 months  
*Storage Temp:* -20°C  
*Storage Conditions:* 10 mM Tris-HCl; 1 mM EDTA @ pH 8.0  
*Specification Version:* PS-N9215S v2.0  
*Effective Date:* 07 Dec 2017

### Assay Name/Specification (minimum release criteria)

**A260/A280 Assay** - The ratio of UV absorption of pCLIPf Vector at 260 and 280 nm is between 1.8 and 2.

**DNA Concentration (A260)** - The concentration of pCLIPf Vector is between 500 and 550 µg/ml as determined by UV absorption at 260 nm.

**Electrophoretic Pattern (Plasmid)** - The banding pattern of pCLIPf Vector on a 1.2% agarose gel is evaluated against a control lot for sharpness and relative intensity as determined by gel electrophoresis using Ethidium Bromide.

**Non-Specific DNase Activity (DNA, 16 hour)** - A 50 µl reaction in 1X NEBuffer 2 containing 2.5 µg of pCLIPf Vector incubated for 16 hours at 37°C results in a DNA pattern free of detectable nuclease degradation as determined by agarose gel electrophoresis.



Date 07 Dec 2017

Derek Robinson  
Director of Quality Control

